


# PRODUCT DATA SHEET

Revision 1	384 Well Cell Culture Plate, White, F-Bottom	
Validfor Item-No.:	1015100/1015104/1015160/1015164/1015102/1015106/ 1015162/1015166	

1.	Description / Specification	
1.1	Description	Cell Culture Plate, 96 well, solid white flat bottom , with vents, alphanumeric well coding, single position lid with condensation rings, sterile 1015100 384 well, Cell culture plate, with white lid, TC treated, sterile 1015104 384 well, Cell culture plate, with white lid, TC treated, sterile 1015160 384 well, Cell culture plate, with clear lid, TC treated, sterile 1015164 384 well, Cell culture plate, with clear lid, TC treated, sterile 1015102 384 well, Cell culture plate, with white lid, TC treated, sterile 1015106 384 well, Cell culture plate, with white lid, TC treated, sterile 1015162 384 well, Cell culture plate, with clear lid, TC treated, sterile 1015166 384 well, Cell culture plate, with clear lid, TC treated, sterile
1.2	Dimensions	Plate : see Customer Drawing 1015100-1015164:Individually Wrapped in Peelable 1015102-1015166:Individually paper-plastic wrapped
1.3	Volume	Max volume: 0.112 ml (mathematical calculated) Working volume: 0.02-0.08 ml Growth area / well: 0.06 cm <sup>2</sup>
1.4	Material / Resin	Plate : PS (Polystyrene) Lid: PS (Polystyrene) The materials for manufacturing are <i>USP Class VI</i>
1.5	Colour	Plate :black Lid: clear/ black
1.6	Sterilization	SAL 10 <sup>-6</sup>
1.7	Quality Control	- Raw Material-Control: physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Expiry date and Lot-No. printed on the outer box

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	1015104/1015164/1015106/ 1015166	1015100/1015160/1015102/ 1015162/
3.1	Pieces / Bag	1	1
3.2	Pieces / Box	50	100
3.3	Lot-No.	DDD YY XXX (day, year, product No.)	
3.4	Other Information	Certificate of Quality to download	

4.	Other Information
4.1	Laboratory use only

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 1 June 2024	Date 1 June 2024	Date 1 June 2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.