PRODUCT DATA SHEET			
Revision 1	150 mm Cell Culture Dish	SAINING LIFE SCIENCE	
Validfor Item-No.:	1025040/1025000/1025140/1025100		

1.	Description / Specification			
1.1	Description	Cell Culture Dish, 150 mm, with vents 1025040 150 mm Cell culture dish, easy-grip type, Non-TC treated, sterile 1025000 150 mm Cell culture dish, Non-TC treated, sterile 1025140 150 mm Cell culture dish, easy-grip, TC treated, sterile 1025100 150 mm Cell culture dish, TC treated, sterile		
1.2	Dimensions	Dish: see Customer Drawing		
1.3	Volume	Max volume: 240 ml Working volume: 30-45 ml Growth area / well: 148 cm²		
1.4	Material / Resin	Plate: PS (Polystyrene) <u>Lid:</u> PS (Polystyrene) The materials for manufacturing are <i>USP Class VI</i>		
1.5	Colour	Dish: clear Lid: clear		
1.6	Sterilization	SAL 10 ⁻⁶		
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification		
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment.		
1.9	Other Information	- For single use only - Expiry date and Lot-No. printed on the outer box		

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	5
3.2	Pieces / Box	12
3.3	Lot-No.	DDD YY XXX (day, year, product No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
4.1	Laboratory use only

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in
Revision 1	Date 1 June 2024	Date 1 June 2024	Date 1 June 2024	this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.