PRODUCT DATA SHEET			
Revision 1	6 Well Cell Culture Plate	SAINING LIFE SCIENCE	
Validfor Item-No.:	1010300/1010304/1010000/1010004/1010002/1010006		

1.	Description / Specification			
1.1 Description		Cell Culture Plate, 6 well, with vents, alphanumeric well coding, single position lid with condensation rings, sterile 1010300 6 well, Cell culture plate, Non-TC treated, sterile 1010304 6 well, Cell culture plate, Non-TC treated, sterile		
		1010000 6 well, Cell culture plate, TC treated, sterile 1010004 6 well, Cell culture plate, TC treated, sterile 1010002 6 well, Cell culture plate, TC treated, sterile 1010006 6 well, Cell culture plate, TC treated, sterile		
1.2	Dimensions	Plate : see Customer Drawing 1010300-1010004:Invidually Wrapped in Peelable 1010002-1010006:Individually paper-plastic wrapped		
1.3	Volume	Max volume: 16.8 ml (mathematical calculated) Working volume: 2 – 5 ml Growth area / well: 9,5 cm ²		
1.4	Material / Resin	Plate : PS (Polystyrene) <u>Lid:</u> PS (Polystyrene) The materials for manufacturing are <i>USP Class VI</i>		
1.5	Colour	Plate : clear Lid: clear		
1.6	Sterilization	SAL 10 ⁻⁶		
1.7	Quality Control	- Raw Material-Control: physical testing - <u>Product-Control:</u> testing of attributive and variable characteristics in accordance with the valid specification		
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment.		
1.9	Other Information	- For single use only - Expiry date and Lot-No. printed on the outer box		

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	1
3.2	Pieces / Box	50/100
3.3	Lot-No.	DDD YY XXX (day, year, product No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
4.1	Laboratory use only

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in
Revision	Date	Date	Date	this document or drawing is confidential and proprietary to Zhejiang
1	07.05.2024	14.05.2024	31.05.2024	Saining Biotechnology Co., Ltd. This
Date	Name	Name	Name	document may not be reproduced for any reason without written permission from
01.06.2024	Chen Anjing	Zhang Heng	Wang Kedong	Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.